

COMAR LIQUID DISPENSING SOLUTIONS

A Full Line of Liquid Dispensing Products, from Oral to Topical: Unparalleled Safety and Accuracy!



55 Years of Quality, Innovation and Service



WHEN LIQUID DISPENSING SAFETY, ACCURACY AND CONVENIENCE REALLY MATTER

Comar's patented liquid dispenser and dosing system is safe and accurate for oral or topical applications. Safety is enhanced through the unique design of its two-piece cylinder and plunger system, without the need for a rubber grommet tip to create the perfect seal. As an "oral syringe", when used with Comar's Press-In-Bottle-Adapter and Bottles, drawing, measuring and dispensing is highly accurate and clean. Accurate dosing with Comar's newly designed and upgraded Dosage Cup is also a unique way to "improve accuracy" in the home or patient care environments!

A Little Bit of Innovation Goes a Long Way

- One-Piece Plunger with Patented Wiper Design
- Custom Calibration Lines and Printing
- Tested, Approved and Used by Major Pharmaceutical Firms
- Accurate, User-Friendly, Frictionless Operation
- Accessories for a Complete Liquid Dispensing System

[YOUR PACKAGING SOLUTION]

THE BENEFITS OF COMAR'S ORAL DISPENSER PRODUCTS

■ INNOVATIVE

- One-Piece Plunger Eliminates Rubber Grommet
- Patented Wiper Design

■ ACCURATE

- Calibration Lines 100% Vision Inspected for Accuracy
- Pressure Testing During Assembly to Ensure Functionality

■ SAFE (One-Piece Plunger Design)

- Elimination of Rubber Grommet and Small Parts Reduces Risk of Choking
- Latex Free – Reduced Extractables
- Tip Prohibits Acceptance of Hypodermic Needles

■ EASY AND COST EFFECTIVE CUSTOMIZATION

- 360° Two-Color Custom Printing (Logos, Lines)
- Custom Colors (Plungers, Translucent Barrels, Tip Caps and Stands)

■ DURABLE AND REUSABLE

- Dishwasher Safe (maintains original print quality and functionality through multiple cycles)
- Maintains Print Quality and Seal Integrity over Extended Shelf Life

The Comar Oral Dispenser has a patented one-piece wiper design and quality standards that are unmatched in the industry.

Customize a Comar Liquid Dispensing Solution with Ease

Comar's Liquid Dispensing Products and Accessories offer convenient and cost effective customization in several ways. With the ability to print 360° we can custom print logos, information and graduation lines. Another great way to help your product stand out is with custom color configurations or accessories. In addition to these easy customization options, we can also help you with custom sizes and configurations.

Comar's experienced Sales and Support Teams are ready to help you create a Customized Liquid Dispensing Product.

Call today to get started: 800-962-6627



Comar Oral Dispenser Products and Accessories

ORAL DISPENSER

- 1, 3, 5, 10 and 15mL* Sizes
- Red, Yellow, Blue, Green and White Plungers
- Clear and Amber, Red, Yellow, Blue, or Green Tinted Barrels

PRESS-IN BOTTLE ADAPTERS

- 20, 24, 28 and 33mm Sizes
- Easy, Spill-Free Prefilling of Oral Dispensers
- Innovative Design Works with Glass or Plastic
- Low Profile Allows use of Standard Closure to Seal Bottle

TIP CAPS AND STANDS

- Seals the Oral Dispenser for Prefilling, Storage and Shipping
- Stand for Convenient One-Hand Dosing

UV AMBER BARRELS

- Meets USP <661> for Light Transmission
- Great for Light Sensitive Liquids, Creams and Gels

*Currently Available from a Unit Tool



Comar's Liquid Dispensing Products Can Be Used for an Extensive Range of Applications.



We also stock a large variety of Liquid Dispensing Products that are readily available.

Visit Comar's Online Catalog at www.comar.com/ecatalog to view Liquid Dispensing Products that are readily stocked. Distributors, please visit www.comar.com/distributors to access the Comar Distributor Resource Center. To learn more about Comar's Liquid Dispensing Products, visit www.comar.com/oraldispensers