



Oral Syringe Systems



Unique Two Piece Oral Syringes.

Full line of liquid oral dispensers and accessories.

The Comar Oral Syringe is a two piece dispenser that uses a patented wiper design instead of a rubber grommet to form a complete seal between the plunger and the barrel. With no small siliconized rubber parts in the tip, the risk of choking and product contamination is greatly reduced. Our oral dispensers are ideally suited for use in dispensing medication or applying topical applications.



Oral Syringe Systems

Tested, approved, and used by major pharmaceutical firms, our Oral Syringe has a patented, one-piece wiper design and meets quality standards that are unmatched in the industry.

Plunger Options

- Material: HDPE
- Standard colors: red, yellow, blue, green, white

Barrel Options

- Material: polypropylene
- Standard colors: clear, amber
- UV amber barrels meet USP <671> for light transmission

Customization

- Logos & lines - 360°, 1 or 2-color custom printing
- Colors - plungers, translucent barrels, tips caps, and stands
- Wrapping - bag/pouch, shrink wrap, wrapped with bottle adapter
- Tip options - for compatibility with other adapters and fittings



Innovative

- One-piece plunger eliminates rubber grommet
- Patented wiper design

Accurate

- User-friendly, frictionless operation
- Calibration lines vision-inspected for accuracy
- Pressure tested during assembly to ensure functionality

Safe

- Reduced choking risk - no rubber grommets or small parts
- Reduced extractables - latex free

Durable & Reusable

- Dishwasher safe
- Maintains print quality and seal integrity over extended usage

Oral Syringes



Item Family	Description	Size (mL)	Tip
41-0129	1/2 mL Oral Syringe	.5	Short
41-0008	1 mL Oral Syringe	1	Short
41-0002	3 mL Oral Syringe	3	Long
41-0001	3 mL Oral Syringe	3	Short
41-0004	5 mL Oral Syringe	5	Long
41-0003	5 mL Oral Syringe	5	Short
41-0006	10 mL Oral Syringe	10	Long
41-0005	10 mL Oral Syringe	10	Short
41-0134	20 mL Oral Syringe	20	Short

Tip Caps, & Stands



Item Family	Description	Resin
14-0209	Oral Syringe Tip Cap	HDPE
14-0878	Oral Syringe Self-Righting Tip Cap	HDPE
14-0223	Oral Syringe No-Choke Tip Cap	HDPE

Bottle Adapters



Item Family	Description	Style	Size (mm)	Resin
22-1230	DoseGuard®, Press-In Valved Bottle Adapter	Press-In	20	LDPE
22-0454	Press-In Bottle Adapter	Press-In	20	LDPE
22-0625	Vented Press-In Bottle Adapter	Press-In Vented	20	LDPE
22-0218	Press-In Bottle Adapter	Press-In	24	LDPE
22-0198	Press-In Bottle Adapter	Press-In	28	LDPE
22-0231	Press-In Bottle Adapter	Press-In	33	LDPE

DoseGuard® Press-In Valved Bottle Adapter

The next level of safety for child-safe liquid medicine packaging and delivery just go better.

- Designed to function in industry standard bottles of HDPE, LDPE, PP and PET
- Adapter securely connects to the container finish and rests on top of the finish, this low profile allows for use with standard closures
- Features a patented valve design that remains closed until the tip of the oral syringe is inserted into the valve
- Valve opening is customizable to specific product viscosities and is easily adaptable to an open hole flow restrictor
- Single material one-piece construction reduces contamination risks

